

Saturation Binding Assay

Charles “Mitch” Boudreaux
Talía Albert

November 6th, 2023

Contents

Overview	1
Determining Protein Concentration	2
Hot Ligand Addition	3
Tables	3
Cold Ligand (Reference) Addition	4
Tables	4
Membrane (Protein) Addition	5
Filtering	5
Formulas	5
Protein Concentration calculation	5
3H Hot Ligand calculation	6
Cold Ligand calculation	6
Membrane Receptor calculation	7
Determining Starting Concentration calculation	7

Overview

Use 6 columns of 96-well shallow plate, can test 2 membrane/protein per plate.
Final Volume in each well is **125 uL**

- 25 uL 3H Hot Ligand
- 25 uL Binding Buffer (BB)/Cold Ligand (Reference)
- 75 uL Membrane (Protein)

Table 1: Plate layout after step 7

	Columns 1/7	Columns 2/8	Columns 3/9	Columns 4/10	Columns 5/11	Columns 6/12
A	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem
B	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem
C	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem
D	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem
E	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem
F	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem
G	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem
H	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL BB 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem	25 uL 3H 25 uL Ref 75 uL Mem

Determining Protein Concentration

1. Resuspend pellet in 12 mL Lysis Buffer (10mM Tris + 5% Sucrose, pH 7.4):
 1. Use 1 mL Lysis Buffer to break up pellet
 2. Fill to ~12 mL with Lysis Buffer
 3. Pool pellets into one tube (if multiple pellets)
2. Perform Bradford Protein Concentration Assay:
 1. Sample Preparation:

- 10 uL Pellet Suspension + 790 uL dH2O + 200 uL Bradford Reagent
- 2. Blank Preparation:
 - 10 uL Lysis Buffer + 790 uL dH2O + 200 uL Bradford Reagent
- 3. Incubate @ RT 10 min
- 4. Measure absorbance @ 595 nm
- 5. Calculate protein concentration (Refer to Formulas section or Spreadsheet)

Hot Ligand Addition

3. Prepare ~15 mL appropriate BB w/BSA (~30 uL) in trough
4. In an empty 96-well shallow plate (Need one column):
 1. Add **330 uL** BB to well H and **165 uL** BB to well A-G
 2. Add 3H-Ligand to Well H (Refer to Formulas section or Spreadsheet)
 3. Perform a Serial Dilution (1:2) of **165 uL** up from well H to A
 4. Remove **25 uL** from Well A for radioactivity counts
 5. Using a multichannel pipettor, Transfer **25 uL** into 6 columns of the Drug Plate

Tables

Table 2: Step 4.1, in separate 96-well plate

	Empty Column
A	165 ul BB
B	165 ul BB
C	165 ul BB
D	165 ul BB
E	165 ul BB
F	165 ul BB
G	165 ul BB
H	330 ul BB

Table 3: Step 4.3, Volume after serial dilution completed

	Empty Column
A	330 ul BB
B	165 ul BB
C	165 ul BB
D	165 ul BB

Empty Column	
E	165 uL BB
F	165 uL BB
G	165 uL BB
H	165 uL BB

Table 4: Plate layout after 3H-Ligand Addition

	Columns 1/7	Columns 2/8	Columns 3/9	Columns 4/10	Columns 5/11	Columns 6/12
A	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H
B	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H
C	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H
D	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H
E	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H
F	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H
G	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H
H	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H

Cold Ligand (Reference) Addition

5. Add **25 uL** BB to first 3 columns
6. Prepare Cold Ligand (Reference):
 1. Usually, add **4 uL** Reference compound into **800 uL** BB into an eppendorf tube
 2. Using a single channel pipettor, add **25 uL** Reference Compound into the last 3 columns

Tables

Table 5: Plate layout after cold ligand (reference) addition

	Columns 1/7	Columns 2/8	Columns 3/9	Columns 4/10	Columns 5/11	Columns 6/12
A	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H	25 uL 3H
	 25 uL BB	25 uL BB	25 uL BB	 25 uL Ref	25 uL Ref	25 uL Ref

	Columns 1/7	Columns 2/8	Columns 3/9	Columns 4/10	Columns 5/11	Columns 6/12
B	25 uL 3H 25 uL BB	25 uL 3H 25 uL BB	25 uL 3H 25 uL BB	25 uL 3H 25 uL Ref	25 uL 3H 25 uL Ref	25 uL 3H 25 uL Ref
C	25 uL 3H 25 uL BB	25 uL 3H 25 uL BB	25 uL 3H 25 uL BB	25 uL 3H 25 uL Ref	25 uL 3H 25 uL Ref	25 uL 3H 25 uL Ref
D	25 uL 3H 25 uL BB	25 uL 3H 25 uL BB	25 uL 3H 25 uL BB	25 uL 3H 25 uL Ref	25 uL 3H 25 uL Ref	25 uL 3H 25 uL Ref
E	25 uL 3H 25 uL BB	25 uL 3H 25 uL BB	25 uL 3H 25 uL BB	25 uL 3H 25 uL Ref	25 uL 3H 25 uL Ref	25 uL 3H 25 uL Ref
F	25 uL 3H 25 uL BB	25 uL 3H 25 uL BB	25 uL 3H 25 uL BB	25 uL 3H 25 uL Ref	25 uL 3H 25 uL Ref	25 uL 3H 25 uL Ref
G	25 uL 3H 25 uL BB	25 uL 3H 25 uL BB	25 uL 3H 25 uL BB	25 uL 3H 25 uL Ref	25 uL 3H 25 uL Ref	25 uL 3H 25 uL Ref
H	25 uL 3H 25 uL BB	25 uL 3H 25 uL BB	25 uL 3H 25 uL BB	25 uL 3H 25 uL Ref	25 uL 3H 25 uL Ref	25 uL 3H 25 uL Ref

Membrane (Protein) Addition

7. Prepare ~4 mL membrane receptor:
 1. Add BB volume to new trough (Refer to Formulas section or Spreadsheet)
 2. Add protein to volume new trough (Refer to Formulas section or Spreadsheet)
 3. Using a multichannel pipettor, add **75 uL** protein dilution into all 6 columns
8. Incubate plates @ RT for 1 hour in drawer

Filtering

Formulas

Protein Concentration calculation

Need:

- OD @ 595 nm

Formula:

$$\text{Protein Concentration (ug/uL)} = \frac{\text{OD@595 nm} - 0.094}{0.503}$$

Where:

- 0.094 and 0.503 were determined from experimental procedure, provided by XP.

3H Hot Ligand calculation

Need:

- Starting Concentration (nM)
- 3H-Ligand Specific Activity (Ci/mmol)

Formula:

$$\text{3H-Ligand Vol (uL)} = \frac{330 \text{ (uL)} * \text{Starting Concentration (nM)} * 5 * 1.2}{\text{Specific Activity (Ci/mmol)}^{-1} * 1000000}$$

Where:

- 330 uL is double the volume of 165 uL, (25 uL * 6 wells * 1.1 overage = 165 uL), so we can perform a serial dilution
- 5 is a Dilution Factor, the final volume in each well is 125 uL, we add 25 uL from the Hot-Ligand plate, (125/25=5)
- 1.2 is a 20% overage
- 1000000 is for unit conversion

Cold Ligand calculation

Usually need **4 uL** Reference compound and **800 uL** BB

Need:

- Concentration of Cold Ligand (Reference) Stock (most are 10 mM)

Formula:

$$\text{Reference Vol (uL)} = \frac{800 \text{ uL} * 10 \text{ uM Final Concentration} * 5}{10000 \text{ uM Starting Concentration}}$$

Where:

- 800 uL is approximate volume we need (8 Wells/Column * 3 Columns * 25 uL/Well * 1.1 overage = 660), use 800 so pulling from reference is easier
- Final concentration of Cold Ligand (Reference) in each well is 10 uM
- 5 is a Dilution Factor, the final volume in each well is 125 uL, we add 25 uL from the Hot-Ligand plate, (125/25=5)
- 10000 uM or 10 mM is starting reference concentration, usually

Membrane Receptor calculation

Need:

- Protein Concentration (ug/uL) (Found from “Determining Protein Concentration”)
- Protein/Well (ug) (XP will provide this)

Formula:

$$\text{Volume of Protein (uL)} = \frac{\text{Protein/Well (ug)} * 4000 \text{ (uL)}}{\text{Protein Concentration (ug/uL)} * 75 \text{ (uL)}}$$

$$\text{Volume of BB (uL)} = 4000 \text{ (uL)} - \text{Volume of Protein (uL)}$$

Where:

- 75 uL is the volume that will be dispensed into each well
- 4000 uL is approximate volume we need (8 Wells/Column * 6 Columns * 75 uL/Well * 1.1 overage = 3960 uL), round to 4000 uL for convenience

Determining Starting Concentration calculation

After obtaining actual radioactive counts (dpm) starting concentration (concentration of 3H-Ligand in Well H) can be determined.

Need:

- Specific Activity (Ci/mmol)
- Actual Counts (dpm)

Formula:

$$\text{Starting Concentration (nM)} = \frac{\text{Actual Counts (dpm)} * 10^9 \text{ (nM/M)} * 2^7}{2.22 * 10^{12} \text{ (dpm/Ci)} * \text{Specific Activity (Ci/mmol)} * 0.125 \text{ (mL/Well)}}$$

Where:

- 10^9 in a unit conversion to nM from M
- 2^7 is for the 7 serial 1:2 dilutions
- $2.22 * 10^{12}$ dpm/Ci is a constant
- 0.125 mL/Well is the final volume in each well